

## Desinfektionsmittelliste für CBC NUCOS

### List of Disinfectants for CBC NUCOS

Validierte Desinfektionsmittel / Validated Disinfectants	Verträgliche Desinfektionsmittel / Compatible Disinfectants
<ul style="list-style-type: none"> <li>Schülke mikrozyd® universal wipes premium</li> </ul>	<ul style="list-style-type: none"> <li>B. Braun Meliseptol® acute Wipes</li> <li>B. Braun Meliseptol® sensitive Wipes</li> <li>Bode Kohrsolin® FF</li> <li>Bode Mikrobac® Tissues</li> <li>Clinell® Peracetic acid wipes</li> <li>Clinell® universal wipes</li> <li>Dr. Schumacher Cleanisept® Wipes Maxi</li> <li>Dr. Schumacher Ultrasol OXY® Wipes</li> <li>Ecolab Incidin™ Active (powder)</li> <li>Ecolab Incidin™ Alcohol Wipe</li> <li>Ecolab Incidin™ OxyWipe S</li> <li>Ecolab Incidin™ Pro</li> <li>Ecolab Incidin™ Rapid</li> <li>Ecolab Sekusept™ aktiv (powder)</li> <li>Hartmann BacilloI® Sensitive Tissues</li> <li>Milton™ Sterilising Tablets</li> <li>Sani-Cloth Disinfectant Wipe</li> <li>Schülke mikrozyd® universal wipes premium</li> <li>Schülke perform® (powder)</li> <li>SoChlor® DST (1.7g NaDCC)</li> </ul>

Unverträgliche Desinfektionsmittel / Incompatible Disinfectants:
<ul style="list-style-type: none"> <li>Bode BacilloI® AF</li> </ul>