



CBC NUCOS

The comfortable silicone tourniquet



Focus Open 2022
Special Mention

+ Comfort for patient and caregiver

Intuitive stasis without finger barriers or pinched skin, thanks to the silicone strap with an innovative honeycomb structure. The ergonomic closure allows one-handed operation, saving caregivers valuable time – a win for comfort and efficiency.

+ Hygienic and durable

The reusable CBC NUCOS replaces single-use products, reducing waste and saving costs. Autoclavable at 134 °C, easy to disinfect with disinfection wipes, and tested for 11,200 uses and 400 sterilization cycles – a long-lasting, hygienic alternative.

+ Advantages

- + Suitable for MRI
- + Classic, sturdy snap-in closure
- + No unintentional opening when under tension
- + Extra button for slow release
- + Durable, high-temperature-resistant material
- + Made in Germany – Internationally patented

Learn more



+ Available in two sizes

CBC NUCOS for adults and CBC NUCOS Kids.




CBC NUCOS Kids: REF 500581
Strap: 12.5 mm x 350 mm



CBC NUCOS: REF 500580
Strap: 25 mm x 500 mm

Easy handling



1 Button for gentle release | 2 Open



1 Button for gentle release | 2 Open

How to use

