



CBC Penlight -
your reliable tool for
precise diagnostics.



Discover the
CBC Penlight!

Simply scan the QR code
or learn more at
www.kimetec.de



CBC Penlight

Diagnostic light
with LED



Kimetec GmbH
Dornierstr. 27
71254 Ditzingen
Germany

Tel.: +49 (0) 7156 / 17602-200
Fax.: +49 (0) 7156 / 17602-500
E-Mail: info@kimetec.de
www.kimetec.de

P2070002EN_2026-02



Efficient, durable and precise

The CBC Penlight from Kimetec ensures precise diagnostics with modern LED technology – ideal for everyday medical use and demanding applications

+ Details and advantages

- + Medical device Class I
- + Approved for eye and throat diagnostics
- + CE marking permanently laser-engraved
- + Classified according to DIN EN 62471 "Exempt Group"
- + Warm white LED – color temperature 3,000 K
- + Light intensity 15,000–31,000 mcd (millicandela)
- + Illumination time approx. 100 hours
- + LED lifetime approx. 100,000 hours
- + Cost-efficient and environmentally friendly
- + Momentary or constant-on switching
- + Replaceable AAA batteries (incl. 2 batteries)
- + Housing made of anodized aluminum – shatterproof and impact-resistant
- + Weight: approx. 50 g
- + Dimensions: Length 142 mm, Ø 11.5 mm
- + Made in Germany



Versatile designs for your style

The CBC Penlight is available in six colors and can be customized with your own logo as a high-quality laser engraving.

black
REF 200201

silver
REF 200250

red
REF 200270

purple
REF 200280

blue
REF 200260

lime
REF 200291



Also available with accessories:

Spatula Holder
REF 2002

